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5. Literatuur

1. ELN - WP10 - Consensual recommendations on preanalytical precautions for the immunophenotyping of leukemia and immunoproliferative disorders
2. Stelzer GT, Marti G, Hurley A, McCoy P Jr, Lovett EJ, Schwartz A. U.S.-Canadian Consensus recommendations on the immunophenotypic analysis of hematologic neoplasia by flow cytometry: standardization and validation of laboratory procedures. *Cytometry* 1997; 30: 214-30.
3. Johansson et al. Guidelines on the use of multicolour flow cytometry in the diagnosis of haematological neoplasms. *Brit J Haematol.* 2014; 165: 455-88.
4. de Jongste AH, Kraan J, van den Broek PD, Brooimans RA, Bromberg JE, van Montfort KA, Smitt PA, Gratama JW. Use of TransFix™ cerebrospinal fluid storage tubes prevents cellular loss and enhances flow cytometric detection of malignant hematological cells after 18 hours of storage. *Cytometry B Clin Cytom.* 2014; 86: 272-9.
5. Renshaw SA, Gupta S, Campos M, Hodes L, Renshaw AA, Gould EW. The addition of RPMI significantly improves the cellularity of cerebrospinal fluid cytology specimens over time. *Cancer Cytopathol.* 2013; 121: 271-4.
6. Bone marrow collection, Youtube
7. Euroflow document referentie <https://euroflow.org/usr/pub/qa.php>
8. NVC verificatie document inklaren antistoffen: NVC en VHL website
9. Keeney M, Chin-Yee I, Weir K, Popma J, Nayar R, Sutherland DR. Single platform flow cytometric absolute CD34+ cell counts based on the ISHAGE guidelines. *International Society of Hematotherapy and Graft Engineering. Cytometry.* 1998; 34: 61-70
10. Sutherland DR1, Keeney M, Illingworth A Practical guidelines for the high-sensitivity detection and monitoring of paroxysmal nocturnal hemoglobinuria clones by flow cytometry. *Cytometry B Clin Cytom.* 2012; 82: 195-208.
11. BAL consensus protocol, NVC werkgroep BAL (www.cytometrie.nl)
12. Brooimans RA, Kraan J, van Putten W, Cornelissen JJ, Löwenberg B, Gratama JW. Flow cytometric differential of leukocyte populations in normal bone marrow: influence of peripheral blood contamination. *Cytometry B Clin Cytom.* 2009; 76: 18-26.
13. <http://www.iqproducts.nl/products/perinatal/fetal-cell-count-kit-diagnosis-of-fetomaternal-hemorrhage/of-fetomaternal-hemorrhage/>
14. Shabnam Tanqri, Validation of Cell-based Fluorescence Assays: Practice Guidelines from the ICSH and ICCS- Part III - Analytical Issues *Cytometry: Part B (Clinical Cytometry)* 84B:291-308 (2013)
15. Bruce H. Davis, Validation of Cell-based Fluorescence Assays: Practice Guidelines from the ICSH and ICCS. - Part II - Preanalytical Issues *Cytometry: Part B (Clinical Cytometry)* 84B:286-290 (2013)
16. Richtlijn paroxysmale nachtelijke hemoglobinurie. Landelijke werkgroep PNH, NvVH, S. Langemeier, Nijmegen 2016